



PRODUCT CATALOG · AMBU

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 1

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

MEDICAL CONSUMABLES

- Neuroline Concentric & Monopolar EMG Needles

Medical Consumables

Neuroline Concentric & Monopolar EMG Needles



Neuroline monopolar needle electrodes

COLOR	ITEM NO.	LENGTH (MM)	LENGTH (IN)	DIAMETER (MM)	GAUGE
Grey	74225-36/40	25	1	0.35	29G
Light green	74238-36/40	38	1.5	0.35	29G
Green	74238-45/40	38	1.5	0.45	26G
Yellow	74250-45/40	50	2	0.45	26G
Red	74275-50/40	75	3	0.50	25G
Grey	74325-36/40*	25	1	0.35	29G
Light green	74338-36/40*	38	1.5	0.35	29G
Green	74338-45/40*	38	1.5	0.45	26G
Yellow	74350-45/40*	50	2	0.45	26G
Red	74375-50/40*	75	3	0.50	25G

Neuroline concentric needle electrodes

COLOR	ITEM NO.	LENGTH (MM)	LENGTH (IN)	DIAMETER (MM)	GAUGE	SHELF LIFE (MO)	RECORDING AREA (MM ²)
White	74025-30/25	25	1	0.30	30G	36	0.02
Grey	74025-45/25	25	1	0.45	26G	36	0.07
Purple	74030-35/25	30	1.2	0.36	28G	36	0.02
Green	74038-45/25	38	1.5	0.45	26G	36	0.07
Yellow	74050-45/25	50	2	0.45	26G	36	0.07
Red	74075-65/25	75	3	0.64	23G	36	0.07

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.