



PRODUCT CATALOG · DEYMED

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 2

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

NEUROLOGY

- DuoMAG TMS / rTMS
- Neurofeedback System

Neurology

DuoMAG TMS / rTMS

DuoMAG MP



Type	Monophasic TMS stimulator (single-pulse)
Pulse energy	High-energy monophasic (~700 J per pulse)
Max repetition rate	Up to 2 pps (0.5 Hz at 100% intensity, single unit)
Mains	100–240 VAC, 50/60 Hz
Sync / control	TTL in/out and/or USB; coil controls; TruScan EEG / TruTrace EMG control
Configurations	MagCart, MagTower, Compact (Compact not compatible with coil cooling)
Variant	DuoMAG MP-EEG for advanced TMS-EEG (research)
Regulatory note	EU MDR certified where registered; investigational device in the USA

DuoMAG XT / XTR DuoMAG XT



Type	Biphasic rTMS / TMS stimulator
Models	XT-10, XT-35, XT-100, XTR
Pulse modes	Single / repetitive (all); burst / TBS on XT-35, XT-100, XTR
Max repetition rate	XT-10: 30 pps; XT-35: 60 pps; XT-100 / XTR: 100 pps
Mains	100–240 VAC, 50/60 Hz
Sync / control	TTL in/out and/or USB; coil controls; TruScan EEG / TruTrace EMG control
Configurations	MagTower, MagCart, Compact (Compact not compatible with coil cooling)
Regulatory note	EU MDR certified where registered; investigational device in the USA

DuoMAG overview

DUOMAG MP	DUOMAG XT / XTR
Pulse shape	Monophasic
Primary mode	Single-pulse
Max rate	Up to 2 pps
Typical use	MEP, high-energy monophasic TMS

DuoMAG XT power classes

XT-10	XT-35	XT-100	XTR
Pulse modes	Single, repetitive	Single, repetitive, burst	Single, repetitive, burst
Max repetition rate	30 pps	60 pps	100 pps

Neurofeedback System



Type	EEG neurofeedback (biofeedback) system
Configurations	FlexiCart (integrated PC + Click N' Go) or Portable (laptop + USB adapter)
Headbox options	2-channel bipolar EEG; 4-channel (2× EEG, 1× EMG, 1× GSR)
Sample rate	256 Hz
Power / isolation	Battery or direct power; optical isolation; inductive charging on cart arm
Connectors	Touch-proof electrode connectors
Software	Dashboard UI; Z-score & custom protocols; FFT; 6 threshold bins; 20+ 3D games
Regulatory note	Not available for sale in the USA; not an EU medical device (no diagnosis/treatment claims in EU)

Configurations & headboxes

FLEXICART	PORTABLE
Form factor	Cart with integrated PC + Click N' Go
Headbox options	2ch or 4ch
2ch headbox	Bipolar EEG (ref / active / ground per channel)
4ch headbox	2× EEG + EMG + GSR

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.