



PRODUCT CATALOG · DM INFORMATION

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 3

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

CARDIOLOGY

- Ambulatory BP Monitor
- Holter Recorder and Analyzer
- Stress Test ECG

Cardiology

Ambulatory BP Monitor



Type	Ambulatory blood pressure monitor (upper arm)
Method	Oscillometry with step deflation
Recording period / memory	1–48 hours; up to 460 complete measurements
BP / pulse range	20–260 mmHg; 30–220 bpm
Accuracy / test duration	±3 mmHg; typically 30–90 seconds per measurement
Display / transfer	LCD (SYS/DIA/pulse + status); USB to PC analysis software
Power	2 × AA (LR06) alkaline; DC 3.0 V
Dimensions / weight	Approx. 116 × 70 × 40 mm; ~252 g including batteries

Holter Recorder and Analyzer



Type	Ambulatory Holter ECG recorder
Leads / channels	3-lead or 12-lead (auto-detected by cable: 4/5/7 vs 10 electrodes)
Recording time	Multi-day ambulatory (typically up to ~5 days on supported configs)
Storage	Removable SD card (≥1 GB class recommended)
Display / controls	LCD screen with navigation and enter keys; realtime ECG channel view
Power	1 × AAA alkaline (new battery per recording)
Data transfer	SD card reader or HDMI-to-USB download cable → CardioScan software
Dimensions / weight	Approx. 88 × 55 × 21 mm; ≤100 g (recorder body)

Stress Test ECG



Type	PC-based exercise stress ECG system
Model	DMI-SE01 (SE02/SE03 variants without treadmill/bike options)
Acquisition	12-lead synchronous sampling; Bluetooth wireless or USB wired modes
Exercise devices	Compatible medical treadmills and ergometers (optional with SE01)
Analysis	ST-segment analysis; dynamic/static/superimposed ECG review; multi-parameter trends
Blood pressure	Manual/automatic stage BP entry; optional stress BP monitor docking
Export / interfaces	PDF, ASCII, SQL, GDT, HL7, DICOM; HIS integration on supported setups

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.