



PRODUCT CATALOG · IMPETO MEDICAL

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 1

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

NEUROLOGY

- SUDOSCAN 2

Neurology

SUDOSCAN 2



Type	Sudomotor / electrochemical skin conductance (ESC) analyzer
Principle	Low-voltage DC stimulus on hand & foot plates; chloride-ion release → ESC
Test duration	Approximately 3 minutes
Patient interface	4 sensor plates (hands + feet); Smart Electrodes or standard plates
Control unit	S2MU touchscreen Windows 10 PC (master unit)
Max patient weight on plates	199 kg / 438 lbs
Operating room temperature	18–35 °C (65–95 °F)
Classification	EU Class IIa (93/42/EEC); FDA Class II (21 CFR 882.1540); type BF applied parts
Standards	IEC 60601-1 Ed. 2 / 3.1; IEC 60601-1-2 Ed. 3 / 4
Key accessories (codes)	APH2, APF2, DSPHD, DSPFD, DSP, PWS2, ASAFE, ELOCK2
Manufacturer	Impeto Medical SAS, Paris, France

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.