



PRODUCT CATALOG · MAGSTIM

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 1

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

NEUROLOGY

- Magstim Rapid Therapy System

Neurology

Magstim Rapid Therapy System

New Rapid 6000-00



Type	rTMS / TMS stimulator (Standard Rapid)
Max frequency	50 Hz
Part numbers	6000-00 (230 V, 3 kVA); 6000-US (115 V, 2.3 kVA)
Interface	15" touchscreen; Windows-based software platform
MEP / EMG	Improved 2-channel MEP with charge-delay artifact reduction
Pulse modes	Biphasic; repetitive; standard and burst (incl. theta burst)
Clearance (US)	Peripheral-nerve stimulation for chronic intractable, post-traumatic, and post-surgical pain (18+)

New Super Rapid 6100-00



Type	rTMS / TMS stimulator (Super Rapid)
Max frequency	100 Hz
Part numbers	6100-00 (230 V, 6 kVA); 6100-US (115 V, 4.6 kVA)
Power architecture	Dual HV power supply unit (PSU)
Interface	15" touchscreen; Windows-based software platform
MEP / EMG	Improved 2-channel MEP with charge-delay artifact reduction
Upgrade	From Standard Rapid (6000-00) via 3014-00 kit + mains-lead filter

New Super Rapid Plus1 6200-00



Type	rTMS / TMS stimulator (Super Rapid Plus1)
Max frequency	100 Hz
Part numbers	6200-00 (230 V, 9 kVA); 6200-US (115 V, 6.9 kVA)
Power architecture	Dual HV PSU + Plus1 combining module
Interface	15" touchscreen; Windows-based software platform
MEP / EMG	Improved 2-channel MEP with charge-delay artifact reduction
Upgrade	From Super Rapid (6100-00) with Plus1 combining module

Model overview

NEW RAPID	NEW SUPER RAPID	NEW SUPER RAPID PLUS1
Part no. (230 V)	6000-00	6100-00
Part no. (115 V)	6000-US	6100-US
Max frequency	50 Hz	100 Hz
Power (230 V)	3 kVA	6 kVA
Power (115 V)	2.3 kVA	4.6 kVA

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.