



PRODUCT CATALOG · MEDICAP

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 1

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

VASCULAR DIAGNOSTICS

- Transcutaneous Oxygen Pressure (PRÉCISE 8008)

Vascular Diagnostics

Transcutaneous Oxygen Pressure (PRÉCISE 8008)



Type	Multi-probe transcutaneous oxygen (tcpO ₂ / TCOM) monitor
Probes	Up to 8 optical fluorescence sensors (simultaneous)
Display	15.1" touchscreen
Measuring range	0 to 2000 mmHg (±10%)
Sensor temperature	Adjustable 40–44°C
Data / connectivity	Internal database; USB export; optional visualization software / report PDF
Dimensions / weight	Approx. 400 × 250 × 170 mm; ~3.9 kg (incl. 8 sensors)
Power	100–230 VAC ±10%; max. 35 VA; protection class II, Type BF, IP22

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.