



PRODUCT CATALOG · MEDICO

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 1

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

MEDICAL CONSUMABLES

- Disposable ECG Electrodes

Medical Consumables

Disposable ECG Electrodes

MLGLT-01 Liquid Gel MLGLT-01



Gel system	Liquid
Backing / adhesive	Foam, enhanced
Dimensions (h x w x t)	57 x 38 x 1 mm
Suggested use	Stress / Holter / resting / monitoring

MSGLT-01G Solid Gel MSGLT-01G



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h x w x t)	57 x 38 x 1 mm
Suggested use	Stress / Holter / resting / monitoring

MSGLT-08G Solid Gel MSGLT-08G



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h x w x t)	45 x 42 x 1 mm
Suggested use	Stress / resting / monitoring

MSGLT-08GRT Solid Gel MSGLT-08GRT



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h × w × t)	45 × 42 × 1 mm
Suggested use	Stress / resting / monitoring

MSGST-37 Pediatric MSGST-37



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h × w × t)	35 × 30 × 1 mm
Suggested use	Resting / monitoring (pediatric)

Stress / Holter foam electrode configurations

PRODUCT CODE / GEL	BACKING / ADHESIVE	DIMENSIONS (H × W × T)	SUGGESTED USE
MLGLT-01, liquid MSGLT-01G, solid*	Foam, enhanced	57 × 38 × 1 mm	Stress / Holter / resting / monitoring
MSGLT-08G, solid	Foam, enhanced	45 × 42 × 1 mm	Stress / resting / monitoring
MSGLT-08GRT, solid	Foam, enhanced	45 × 42 × 1 mm	Stress / resting / monitoring
MSGST-37, solid†	Foam, enhanced	35 × 30 × 1 mm	Resting / monitoring

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.