



PRODUCT CATALOG · HYPERBARIC THERAPY

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 5

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

OXYHEALTH

- [Fortius 420 Hyperbaric Chamber](#) Hyperbaric Therapy
- [Quamvis 320 Hyperbaric Chamber](#) Hyperbaric Therapy
- [Respiro 270 Hyperbaric Chamber](#) Hyperbaric Therapy
- [Solace 210 Hyperbaric Chamber](#) Hyperbaric Therapy
- [Vitaeris 320 Hyperbaric Chamber](#) Hyperbaric Therapy

OxyHealth

HYPERBARIC THERAPY

Fortius 420 Hyperbaric Chamber



Type	Hard-shell hyperbaric chamber
Operating pressure	Up to 3.0 ATA / 29 psi (US); EU/UK wellness 1.3–2.0 ATA typical
Internal diameter	106.68 cm (42 in)
Length	231.14 cm (91 in)
Door diameter	66.80 cm (26.3 in)
Chamber weight	567 kg (1,250 lbs)
Compressor	Double-headed oil-less 1/3 hp; 0.01 µm filtration; ~16.8 kg
Occupancy	Comfortably fits two adults

HYPERBARIC THERAPY

Quamvis 320 Hyperbaric Chamber



Type	Soft-shell mild hyperbaric chamber
Operating pressure	1.3 ATA (4 psi)
Classification (EU & UK)	Wellness non-medical device
Inflated diameter	81.28 cm (32 in)
Length	254 cm (100 in)
Bladder / frame	17.23 kg (38 lbs) / 15.42 kg (34 lbs)
Mattress / bolsters	5.89 kg (13 lbs) / 3.17 kg (7 lbs)
Compressor	Double-headed oil-less 1/3 hp; 0.01 µm filtration; ~16.8 kg

HYPERBARIC THERAPY

Respiro 270 Hyperbaric Chamber



Type	Soft-shell mild hyperbaric chamber
Operating pressure	1.3 ATA (4 psi)
Classification (EU & UK)	Wellness non-medical device
Inflated diameter	68.58 cm (27 in)
Length	223.52 cm (88 in)
Bladder / frame / mattress	12.24 kg / 7.25 kg / 4.53 kg
Bolsters	3.17 kg (7 lbs)
Compressor	Double-headed oil-less 1/4 hp; 0.01 µm filtration; ~10.9 kg

HYPERBARIC THERAPY

Solace 210 Hyperbaric Chamber



Type	Soft-shell mild hyperbaric chamber
Operating pressure	1.3 ATA (4 psi)
Classification (EU & UK)	Wellness non-medical device
Inflated diameter	53.34 cm (21 in)
Length	233.68 cm (92 in)
Bladder / frame / mattress	10.88 kg / 5.44 kg / 2.72 kg
Bolsters	3.17 kg (7 lbs)
Compressor	Double-headed oil-less; 0.01 µm filtration; ~10.9 kg

HYPERBARIC THERAPY

Vitaeris 320 Hyperbaric Chamber



Type	Soft-shell mild hyperbaric chamber
Operating pressure	1.3 ATA (4 psi)
Classification (EU & UK)	Wellness non-medical device
Inflated diameter	81.28 cm (32 in)
Length	233.68 cm (92 in)
Bladder weight	16 kg (35 lbs)
Frame / mattress / bolsters	7 kg / 4 kg / 3 kg
Compressor	Double-headed oil-less 1/3 hp; 0.01 µm filtration; ~17 kg

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.