



PRODUCT CATALOG · MEDICAL CONSUMABLES

ICare Life Medical Devices, Inc.

Medical devices for hospitals and clinics in the Philippines

Generated August 8, 2026

Active product lines 9

Website icarelifemedical.com

Address Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email sales@icarelifemedical.com

Phone +63 (02) 8741 8043 · +63 917 854 9704

Office hours Monday–Friday • 8:00 AM – 5:00 PM

Contents

NIHON KOHDEN

- [EEG Disk Electrodes \(Collodion\)](#) Medical Consumables
- [Elefix V EEG Paste](#) Medical Consumables
- [Skin Pure Skin Preparation Gel](#) Medical Consumables

AMBU

- [Neuroline Concentric & Monopolar EMG Needles](#) Medical Consumables

HONGKONG BEST MEDICAL

- [Patient Monitoring & Diagnostic Accessories](#) Medical Consumables

MEDICO

- [Disposable ECG Electrodes](#) Medical Consumables

TELE-PAPER

- [ECG & Medical Chart Paper](#) Medical Consumables
- [Stress-Test Chart Paper](#) Medical Consumables
- [Ultrasound / Video Printer Paper](#) Medical Consumables

Nihon Kohden

MEDICAL CONSUMABLES

EEG Disk Electrodes (Collodion)



Manufacturer code	NE-134A
Disk diameter	10 mm silver disk
Lead length	60 in (approx. 1.5 m)
Connector	1.5 mm protected connectors
Packaging	12 electrodes per pack

MEDICAL CONSUMABLES

Elefix V EEG Paste



Elefix V pack sizes

PACK	MANUFACTURER CODE	SCOPE
400 g × 3 jars	ZV-401E03	EEG
180 g × 10 tubes	ZV-181E10	EEG
180 g × 2 tubes	ZV-181E02	EEG

MEDICAL CONSUMABLES

Skin Pure Skin Preparation Gel



Volume	135 g × 2 tubes
Description	Skin preparation gel for stable electrode contact
Scope of application	ECG and EEG
Manufacturer code	YZ-0019
Order code	F020

Ambu

MEDICAL CONSUMABLES

Neuroline Concentric & Monopolar EMG Needles



Neuroline monopolar needle electrodes

COLOR	ITEM NO.	LENGTH (MM)	LENGTH (IN)	DIAMETER (MM)	GAUGE
Grey	74225-36/40	25	1	0.35	29G
Light green	74238-36/40	38	1.5	0.35	29G
Green	74238-45/40	38	1.5	0.45	26G
Yellow	74250-45/40	50	2	0.45	26G
Red	74275-50/40	75	3	0.50	25G
Grey	74325-36/40*	25	1	0.35	29G
Light green	74338-36/40*	38	1.5	0.35	29G
Green	74338-45/40*	38	1.5	0.45	26G
Yellow	74350-45/40*	50	2	0.45	26G
Red	74375-50/40*	75	3	0.50	25G

Neuroline concentric needle electrodes

COLOR	ITEM NO.	LENGTH (MM)	LENGTH (IN)	DIAMETER (MM)	GAUGE	SHELF LIFE (MO)	RECORDING AREA (MM ²)
White	74025-30/25	25	1	0.30	30G	36	0.02
Grey	74025-45/25	25	1	0.45	26G	36	0.07
Purple	74030-35/25	30	1.2	0.36	28G	36	0.02
Green	74038-45/25	38	1.5	0.45	26G	36	0.07
Yellow	74050-45/25	50	2	0.45	26G	36	0.07
Red	74075-65/25	75	3	0.64	23G	36	0.07

HongKong Best Medical

MEDICAL CONSUMABLES

Patient Monitoring & Diagnostic Accessories



Medico

MEDICAL CONSUMABLES

Disposable ECG Electrodes

MLGLT-01 Liquid Gel MLGLT-01



Gel system	Liquid
Backing / adhesive	Foam, enhanced
Dimensions (h x w x t)	57 x 38 x 1 mm
Suggested use	Stress / Holter / resting / monitoring

MSGLT-01G Solid Gel MSGLT-01G



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h x w x t)	57 x 38 x 1 mm
Suggested use	Stress / Holter / resting / monitoring

MSGLT-08G Solid Gel MSGLT-08G



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h x w x t)	45 x 42 x 1 mm
Suggested use	Stress / resting / monitoring

MSGLT-08GRT Solid Gel MSGLT-08GRT



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h × w × t)	45 × 42 × 1 mm
Suggested use	Stress / resting / monitoring

MSGST-37 Pediatric MSGST-37



Gel system	Solid
Backing / adhesive	Foam, enhanced
Dimensions (h × w × t)	35 × 30 × 1 mm
Suggested use	Resting / monitoring (pediatric)

Stress / Holter foam electrode configurations

PRODUCT CODE / GEL	BACKING / ADHESIVE	DIMENSIONS (H × W × T)	SUGGESTED USE
MLGLT-01, liquid MSGLT-01G, solid*	Foam, enhanced	57 × 38 × 1 mm	Stress / Holter / resting / monitoring
MSGLT-08G, solid	Foam, enhanced	45 × 42 × 1 mm	Stress / resting / monitoring
MSGLT-08GRT, solid	Foam, enhanced	45 × 42 × 1 mm	Stress / resting / monitoring
MSGST-37, solid†	Foam, enhanced	35 × 30 × 1 mm	Resting / monitoring

Tele-Paper

MEDICAL CONSUMABLES

ECG & Medical Chart Paper



ECG & chart paper sizes

DESCRIPTION	SIZE (MM)	USE
ECG paper SM5030/16R3	50 × 30	ECG roll
ECG paper SM5020/16R3	50 × 20	ECG roll
thermal paper SM6330/16R3	63 × 30	ECG roll
Hellige Marquette HM210295/150RS	210 × 295	ECG Z-fold
Nihon Kohden ECG thermal paper NK210140/142RS	210 × 140 × 142 sh	ECG Z-fold
Nihon Kohden ECG thermal paper NK210140/200RS	210 × 140 × 200 sh	ECG Z-fold
Nihon Kohden Evertrace heat-sensitive paper	63 × 30	ECG roll
ECG roll CSO SM8025/16R3	80 × 25	ECG roll
ECG roll CSO SM5030/16BS3	50 × 30	ECG roll
Kenz R112-27Z/272-C2 (K11290/300RS)	112 × 90	ECG Z-fold
Recording chart paper Comen STAR5000E	112 × 90	CTG / fetal
Edan/Corometrics AAO4305/CAO4305 (CO15290/150R)	152 × 90	CTG / fetal

MEDICAL CONSUMABLES

Stress-Test Chart Paper



Stress-test paper sizes

DESCRIPTION	SIZE (MM)	USE
Marquette thermal paper MQ215280/300R	215 × 280	Stress test
Nihon Kohden ECG thermal paper (sheet format)	210 × 295 × 250 sh	Stress test / ECG

MEDICAL CONSUMABLES

Ultrasound / Video Printer Paper



Ultrasound / video printer paper sizes

DESCRIPTION	SIZE (MM)	USE
Video printer paper SM11018 (Type V) glossy	110 × 18	Ultrasound
Video printer paper SM110HD	110 × 20	Ultrasound
Hellige Marquette HM8090/280RS	80 × 90	Ultrasound
Video printer paper SM110S	—	Ultrasound

Contact

ICare Life Medical Devices, Inc.

Zen Jose Bldg., 58 Ninong Isidro St., Brgy. 169, Caloocan City 1422

Email: sales@icarelifemedical.com

Phone: +63 (02) 8741 8043 / +63 (02) 7718 4916 / +63 917 854 9704

Web: <https://icarelifemedical.com>

ICare Life Medical Devices, Inc. is an authorized importer and distributor of medical devices and healthcare technologies.
Availability and specifications may change without notice.